



4200 SERIES PNEUMATIC STOOLS

PIONEERING OPHTHALMIC TECHNOLOGY

4200

SERIES PNEUMATIC STOOLS

Finished for Durability

- Select from 14 easy-to-maintain vinyl colours. Other vinyl colours available by request.
- All exposed metal, with exception of the aluminum base, is chrome-plated steel.
- Base and legs finished in tough, black epoxy enamel.
- Hooded dual-wheel casters.
- 5 legs for stability.

Optional Foot Ring

Order a foot ring with your stool, or at a later date.

Tech Specs

Base diameter: 22.5"



| | 4240 | 4246 |
|-------------------------------|------------------|------------------|
| Seat adjustment | Pneumatic Assist | Pneumatic Assist |
| Seat diameter | 38 cm | 38 cm |
| Seat construction | 5 cm polyfoam | 5 cm polyfoam |
| Seat swivels | 360° | 360° |
| Five legs for extra stability | Yes | Yes |
| Height range | 49.5 - 68.6 cm | 49.5 - 68.6 cm |
| Shipping weight | 10.4 kg | 11.8 kg |