



Tradition and Innovation - Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability. Anticipating trends to improve the quality of life.

7000

DUAL-PURPOSE EXAM CHAIR AND MINOR PROCEDURE TABLE

7000

DUAL-PURPOSE EXAM CHAIR AND MINOR PROCEDURE TABLE



Ergonomic Features

Dual-side operation allows you to control chair from both sides while maintaining patient contact.

Adjustable concave headrest provides comfortable head and neck support, and locks with a single hand.

Upholstered armrests lock in place and release on both chair sides for security and comfort.

Large hinged footrest supports patient's weight and tilts out of the way to save space.

Corded footswitch provides hands-free operation of chair base and top. Dual rotation locks allow you to swivel chair 330°.

Designed for Minor Procedures

The 7000 chair reclines nearly flat for minor procedures. Non-articulating pillow headrest and articulating surgical headrest are optional.

Power Hydraulic base

Square pedestal base provides crucial patient stability and easy accessibility. Finger switches and a corded footswitch raise, lower or recline chair into any position without breaking patient contact. High and low base models available.

Ultra-Durable Membrane Switches

Programmable membrane switches generally last longer than mechanical switches, whilst being backlit is optimal for darkened conditions. By running at low voltages, these switches keep patient safety in mind.

Finishes and Upholstery

- Hand-sewn upholstery covers high-density foam designed to last for years.
- Removable cushions for easy cleaning or replacement.
- Designer-style cushions standard.
- Base available in #20 Pearl.
- Choose from 14 standard vinyl colours and many custom selections.

Dimensions and Specifications

Base	571.5 cm x 61 cm
Overall width	61 cm
Overall depth	91.5 cm (excluding headrest)
Table position	152.4 cm
Seat	48.3 cm wide by 47 cm deep
Armrest distance	48.3 cm

CHOICE OF TWO BASE HEIGHTS:

High	56 cm lowered, 846 cm raised
Low	51 cm lowered, 71 cm raised
Shipping weight	148 kg
Standard electrical	115 volts, 60/50 Hz, 6.3 amperes total
	230 volts, 50 Hz, 4.0 amperes (Optional)

